## TANITA

# MC 780 MA S

### Segmental Multi Frequency Body Composition Analyser

#### **Key Features**

- Full segmental body composition measurements performed in under 20 seconds using the latest multi frequency BIA technology.
- Interactive stand-alone unit allows clients to quickly take a measurement without specialist assistance.
- Extra large LED dual display show analysis in an easy-to-read format.
- Results automatically stored on the SD Card, sent to a PC or transferred to a Pictbridge printer to generate a full consultation sheet for further discussion.
- Data can be sent wirelessly to a PC or Windows Tablet using a Bluetooth adaptor.
- GMon Pro Software compatible, allowing trend analysis, health risk assessments and full data management.
- Max weight capacity 270kg with 100g accuracy
- Accuracy grade: NAWI Class III, MDD Class II-a

#### Measurements

- Body fat %
- Fat mass kg,
- Fat free mass kg
- Muscle mass kg
- Total Body Water %
- Extra Cellular Water kg
- Intra Cellular Water kg
- Body mass index
- Bone mass kg
- Physique rating
- Visceral fat rating
- ECW/TBW ratio
- Basal Metabolic Rate kcal
- Phase Angle
- Basal Metabolic Rate indicator
- Metabolic Age







- - www.tanita.eu

each leg, arm and trunk:
Fat percentage
Fat distribution rating

**O** Segmental readings for

- Muscle mass rating
- Muscle mass kg
- Muscle mass balance
- Leg Muscle Score
- Body Fat Distribution and Reactance/resistance